

Release Notes MatchSource – January 29, 2024

Feature/Enhancement	Change/Value Summary
<p>User-driven Incident Creation</p>	<ul style="list-style-type: none"> • When Technical Difficulties page/pop-up is displayed you can click on 'Submit Incident' to send a ticket directly from MatchSource to the NMDP Service Desk. • Details included in incident submission will include: • Notes section is included for additional context. • Feature allows for faster submission of incidents and includes details necessary for quick resolution from NDMP Service Desk.
<p>Negative HLA Antibody Profile</p>	<ul style="list-style-type: none"> • You can upload a Negative Antibody Report on the 'HLA Antibody Specificity Profile' tab • You can view Negative Antibody Report names in the Virtual Crossmatch (VXM) Filter • You cannot Highlight or Filter on the negative profile • Enhancement allows for viewing all HLA Antibody reports in one location, simplifying transplant source selection
<p>Workup Form: Cryopreservation Shipment</p>	<ul style="list-style-type: none"> • When you plan to use Cryopreservation for either Marrow or PBSC Collection at an AC/CC or NMDP location AND you choose to send a second set of samples to a different location <ul style="list-style-type: none"> • Day of collection samples will ship directly to the TC. If TC cannot receive samples separately from the product, the samples will ship at the same temperature as the product inside the dry shipper container. • This information will display on: <ul style="list-style-type: none"> • PDF exports of a Workup Request form • Order Fulfillment for Workup Request form <ul style="list-style-type: none"> • Original Workup Request tab • Working Workup Request tab • Activated Workup Request tab • Enhancement provides information about cryopreservation shipment/transport limitations earlier in the Workup process to prevent unnecessary back-and-forth communication and rework. • Enhancement aligns with FDA and Customs guidelines for cryopreservation transportation/shipment.
<p>Workup Form: Donor Plasma selection</p>	<ul style="list-style-type: none"> • You will no longer see the 'Donor Plasma' sub-section in the Marrow Collection section of the Workup Form • This information will no longer display on: <ul style="list-style-type: none"> • Marrow Collection section on an exported PDF (in progress and submitted) • Marrow Collection section on the Order Fulfillment UI <ul style="list-style-type: none"> • Original Workup Request tab • Working Workup Request tab

	<ul style="list-style-type: none">• Activated Workup Request tab• Enhancement provides better alignment with network collection center capabilities.
Unborn Patient Weight	<ul style="list-style-type: none">• You are not able to enter the following characters as weight for for an unborn child:<ul style="list-style-type: none">○ 0○ Less than 1 kg○ Decimal value○ Alpha or special characters• The patient weight field can be left blank.• Enhancement allows for more accurate reflection of unknown patient weight for unborn patient.